



Ultrasonic Cleaner Comparison

Doc V: 2.8



Model	Controls	Part Number	Basket Part	Tank Dimensions (mm)			Internal Bath Dimensions (mm)			Basket Dimensions (mm)			Rec'd Fill Level (L)	Ultrasonic Freq (KHz)	Ultrasonic Power (W)	Heater Y / N	Heating (W)	Drain Y/N	Transducers	Mains Power Supply
				L	W	H	L	W	H	L	W	H								
Budget Range																				
0.7 Litre	Switch	US-CU-BU-0.7L	-	-	-	-	140	73	63	146.5	81.5	50	0.4	40	50	N	-	N	1	220V AC
1.3 Litre	Digital	US-CU-4800	Included	-	-	-	160	125	65	-	-	-	-	40	60	N	-	N	1	220V AC
Stainless Steel Range																				
1.3 Litre	Digital	US-CU-DI-1.3L	US-BA-1.3L	180	110	180	140	127	63	140	118	50	0.7	40	50	Y	100	N	1	220V AC
2 Litre	Switch	US-CU-BU-2L	US-BA-2L	180	170	240	138	124	100	140	118	75	1.4	40	50	Y	100	N	1	220V AC
2 Litre	Digital	US-CU-DI-2L	US-BA-2L	180	170	240	138	124	100	140	118	75	1.4	40	50	Y	100	N	1	220V AC
2 Litre	Digital/Cavitek	US-CU-CA-2L	US-BA-2L	180	170	240	138	124	100	140	118	75	1.4	40	50	Y	100	N	1	220V AC
3 Litre	Dial	US-CU-BU-3L	US-BA-3L	265	170	240	238	124	100	235	133	80	2.3	40	100	Y	100	N	2	220V AC
3 Litre	Digital	US-CU-DI-3L	US-BA-3L	265	170	240	238	124	100	235	133	80	2.3	40	100	Y	100	N	2	220V AC
3 Litre	Digital/Cavitek	US-CU-CA-3L	US-BA-3L	265	170	240	238	124	100	235	133	80	2.3	40	100	Y	100	N	2	220V AC
6 Litre	Dial	US-CU-BU-6L	US-BA-6L	370	180	315	275	130	150	270	120	115	5.2	40	150	Y	300	Y	3	220V AC
6 Litre	Digital/Cavitek	US-CU-CA-6L	US-BA-6L	370	180	315	275	130	150	270	120	115	5.2	40	150	Y	300	Y	3	220V AC
9 Litre	Dial	US-CU-BU-9L	US-BA-9L	370	270	315	280	220	150	280	230	115	8.5	40	200	Y	300	Y	4	220V AC
9 Litre	Digital	US-CU-DI-9L	US-BA-9L	370	270	315	280	220	150	280	230	115	8.5	40	200	Y	300	Y	4	220V AC
9 Litre	Digital/Cavitek	US-CU-CA-9L	US-BA-9L	370	270	315	280	220	150	280	230	115	8.5	40	200	Y	300	Y	4	220V AC
10 Litre	Dial	US-CU-BU-10L	US-BA-10L	550	165	310	505	135	150	470	110	100	7	40	200	Y	400	Y	4	220V AC
10 Litre	Digital/Cavitek	US-CU-CA-10L	US-BA-10L	550	165	310	505	135	150	470	110	100	7	40	200	Y	400	Y	4	220V AC
13 Litre	Digital/Cavitek	US-CU-CA-13L	US-BA-13L	-	-	-	330	300	150	315	290	115	11.4	40	300	Y	400	Y	6	220V AC
20 Litre	Digital/Cavitek	US-CU-CA-20L	US-BA-20L	-	-	-	495	295	150	465	270	110	17	40	400	Y	500	Y	8	220V AC
27 Litre	Dial	US-CU-BU-27L	US-BA-27L	530	330	385	495	295	195	465	270	125	23	40	500	Y	500	Y	10	220V AC
27 Litre	Digital	US-CU-DI-27L	US-BA-27L	530	330	385	495	295	195	465	270	125	23	40	500	Y	500	Y	10	220V AC
27 Litre	Digital/Cavitek	US-CU-CA-27L	US-BA-27L	530	330	385	495	295	195	465	270	125	23	40	500	Y	500	Y	10	220V AC
XL Range																				
36 Litre	Dial	US-CU-XL-36L	Included	570	455	360	510	300	250	450	280	140	-	28/40	720	Y	800	Y	14	220V AC
57 Litre	Dial / Digital	US-CU-XL-57L	Included	700	550	670	500	350	360	450	310	170	-	28/40	1.08K	Y	3K	Y	18	220V + 3 Phase
90 Litre	Dial / Digital	US-CU-XL-90L	Included	790	610	710	600	400	400	555	355	210	-	28/40	1.44K	Y	3K	Y	24	220V + 3 Phase
99 Litre	Dial	US-CU-XL-99L	Included	930	380	920	800	250	620	775	220	440	-	40	1.44K	Y	3K	Y	24	220V + 3 Phase
112 Litre	Dial / Digital	US-CU-XL-112L	Included	850	650	710	650	450	400	620	410	210	-	28/40	1.8K	Y	4K	Y	30	220V + 3 Phase
144 Litre	Dial	US-CU-XL-144L	Included	1140	600	670	1000	400	360	960	355	200	-	28/40	1.8K	Y	5K	Y	30	220V + 3 Phase
145 Litre	Dial / Digital	US-CU-XL-145L	Included	900	700	770	700	500	445	680	480	310	-	28/40	2.16K	Y	6K	Y	44	220V + 3 Phase
206 Litre	Dial / Digital	US-CU-XL-206L	Included	960	740	810	750	550	510	720	520	260	-	28/40	2.52K	Y	6K	Y	50	220V + 3 Phase
288 Litre	Dial / Digital	US-CU-XL-288L	Included	1070	920	860	800	600	600	760	560	410	-	28/40	3.6K	Y	6K	Y	72	2x 220V + 3 Phase

To view the full range of tanks and cleaning solutions please visit: www.allendale-ultrasonics.co.uk

